

Frequently asked questions about nasopharyngeal air

What is a nasal aspirator, what effect does it have, what types, what do you notice when using it? Let's find out immediately about this device to help treat common respiratory diseases!

In this article, TipsMake.com will answer common questions about nasal aspirator so you can better understand the product very well in the process of treating these respiratory diseases!

Using a **nasal aspirator** or **nasopharyngeal** nebulizer is an effective alternative to drinking, injecting conventional medicine to treat sore throats, rhinitis, sinusitis . due to dust and pollution. environment, allergens such as pollen, animal hair, weather changes . Find out what you need to know about this device to choose the right product and use the product properly !

Things to know about nose and throat gas machine

1. What is nasopharyngeal nebuliser?
2. Why should a nebulizer be used?
3. Classifying nebulizer
4. What should be noted when buying a nasal aspirator?
5. What should be noted when using a nebulizer?
6. What are the brand of nasopharyngeal machines?
7. How much does the nasopharyngeal gas cost?
8. Buy genuine nasal aspirator, where is it cheap?



The nasopharyngeal air device is suitable for children.

What is nasopharyngeal nebuliser?

A nasal nasal aspirator is a device that helps deliver drugs to respiratory organs in the form of **tiny mist** , which acts directly on infected parts to treat or support the treatment of acute respiratory diseases. and calmness such as *rhinitis, sore throat, laryngitis, trachitis, bronchitis, sinusitis* .

Why should a nebulizer be used?

Should the use of nasopharyngeal aerosol is not a question that people with respiratory diseases, families with young children (who are susceptible to respiratory diseases when changing weather) are interested. So let's see what effect does the nasopharynx have!

Easy to absorb drugs, work quickly, shorten treatment time

When using nebulizer, the drug is converted into tiny mist and then through the nasopharynx, the drug is taken directly to the respiratory organs that need to be affected by inflammation. This fine mist form helps the medicine be introduced into the body faster, can bring the drug directly to the places where other methods are not accessible. Instead of waiting for the medicine to have time to absorb and the drug can be affected by gastric and intestinal fluid . as impact when taking, injecting, aerosol helps the medicine almost immediately work when the user inhales the mist in, thus greatly shortening the duration of treatment.

In addition, besides smoking, users can also use a nasal nasal aspirator with physiological saline to clean the nose and throat, prevent respiratory diseases and support the process of better treatment.



A nasal aspirator helps to treat more respiratory diseases.

Limit side effects

Compared with conventional oral or injectable treatment, aerosol nasal spray will help the drug be taken directly to the diseased tissue, less absorbed into the blood, so it will avoid drugs acting on the whole body, limiting many effects. Side effects of the drug. Especially when using a sore throat instead of taking a direct medicine will limit the side effects causing stomach ulcers, osteoporosis .

Save on treatment costs

When using a nasopharyngeal device, because the drug is absorbed completely through the inhalation route, it will avoid wasting drugs. At the same time, using a home appliance also saves treatment costs more than going to clinics so that the nasal spray can often cost up to 100,000 VND each time.

Easy to use, suitable for many subjects

The nasopharyngeal instrument has a simple structure, easy to use, users can manipulate at home quickly and conveniently. It is suitable for both elderly and children.



Nebulizer has a simple structure, easy to use.

Classifying nebulizer

Currently there are 3 types of aerosol in the market:

1. **Machine used compressed air to create dew seeds:** can be used with suspension, less damaged, durable, affordable. However, this type of device is noisy when operating, providing slower drugs and when the patient breathes, a part of the drug can escape.
2. **Ultrasound machine to create dew seeds:** can provide drugs faster, less noise so patients feel comfortable during the treatment process, and when the patient breathes, the aerosol particles are not escaped. However, this machine is not suitable for drugs with high viscosity or suspension and higher cost than compressed air.
3. **The machine uses ultrasonic vibrating film to create dew seeds:** providing the drug quickly, the treatment effect is much better than the above 2 types, smooth operation, no noise. However, during use, it is necessary to clean the filter net regularly, carefully and the cost of equipment is high, so it is not suitable for use in the home, often used in clinics and hospitals.



Omron NE-C801 nasal aspirator is a nasal aspirator.

What should be noted when buying a nasal aspirator?

Since aerosol generators are a medical and healthcare device, users need to buy genuine products, with clear labels and origins to ensure quality machines, which are effective in the process of conditioning. treatment and does not cause adverse health effects.

Users also need to choose the nasal aspirators with:

1. The fine mist size allows the drug to penetrate into the alveoli, providing the best treatment effect.
2. Has great durability, long life for long-term use.
3. Compact size, convenient to carry people or carry away.

What should be noted when using a nebulizer?

In the process of using the nasal aspirator, users should note some points as follows:

1. Need to see a specialist for examination and get advice on drugs, dosage, smoking time to bring the best treatment effect. Note never pour the medicine over the maximum line.
2. Patients need to keep their posture upright and comfortable during the smoking process. If in the process of suffering from asthma attacks, coughing . then stop immediately to avoid unintended injuries and complications.
3. Regularly clean the device to remove dirt, bacteria, mold . and periodically replace the air duct system to ensure effective nasal spray. Absolutely do not clean the machine by boiling the air duct system, the mask or any part of the machine.

What are the brand of nasopharyngeal machines?

On the market today there are quite a number of brands producing aerosol machines, including the leading brands from Japan and European countries such as Omron nasal aspirator (Japan), gas machine Beurer, air-conditioner Owgels (Germany), Philips throat machine (Netherlands), Laica nasopharyngeal air (Italy) . with advanced and modern production lines, are consumers in difficult markets Information used.



Beurer IH21 nasal aspirator is manufactured in Germany.

How much does the nasopharyngeal gas cost?

Two types of compressed air and ultrasonic devices are commonly used in families with a price ranging from 520,000 to 1,900,000 VND .

Users should consider to buy products with the most suitable price, avoid buying too cheap products that are usually floating, poor quality or counterfeit goods.

Hopefully the above information will be helpful to you in the process of buying a nasal aspirator for your family to help prevent and treat respiratory diseases effectively.

So buy a genuine nasal aspirator, where is it cheap?

META.vn is a unit with over 10 years of experience in providing medical equipment and health care products, committed to genuine sales, warranty and the best prices in the market.

If you need to buy a nasal spray, you can order online on the META.vn website or call **Hanoi** hotline **024.35.68.69.69 - TP. HCM 028.38.33.33.66** for free advice and support for the fastest purchase.

More reference : Top 5 nasal nasal aspirators today

You finished reading the article "**Frequently asked questions about nasopharyngeal air**" edited by the [TipsMake](#) team. We hope this article has provided you with many useful tech tips and tricks. You can search for similar articles on tips and guides. Thank you for reading and for following us regularly.